

the PicoLazer

\$59,995

PICOLAZER ND:YAG LASER
1064nm & 532nm

- Can be used to treat tattoos, benign pigmented lesions & skin revitalization
- Offers high performance speed & comfort
- PicoLazer energy is "non-thermal" and is delivered using a photo-mechanical effect with high-peak powers and the ultra-short pulse durations that shatter pigment without causing thermal damage.
- 450-Pico second pulse duration



ROHRER AESTHETICS

Call or Visit Us Online Now **205.940.2200** www.RohrerAesthetics.com

SPECIFICATIONS

System Type	1064nm / 532nm
Wavelength	1064nm, 532nm
Output	500mJ(1064nm), 250mJ(532nm)
Pulse Duration	450ps
Peak Power	1.33GW(1064nm), 0.8GW(532nm)
Spot Size	3 to 10mm
Repetition Rate	Single, Up to 8Hz
Handpiece	Collimation H/P
Dimensions	25"(W) x 26"(D) x 46"(H)
Weight	176 lbs



ROHRER AESTHETICS

Call or Visit Us Online Now **205.940.2200** www.RohrerAesthetics.com

